
Goat Anti-Mouse IgG (H+L), HRP

Catalog Number: YFSA01; Size: 100 μ L

Background:

Yfxbio® secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA.

We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application.

Product Information:

Product name: Goat Anti-Mouse IgG (H+L), HRP

Product Categories: Secondary Antibody

Form: Liquid

Tested applications: WB, IHC-P, ELISA

Dilution Ratio: WB: 1:5000-1:20000; IHC: 1:500-1:1000; ELISA: 1:5000-1:20000

Species reactivity: Mouse

Host Species: Goat

Conjugation: HRP

Purity: Affinity purification

Immunogen: Mouse IgG

Storage:

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Services & Support:

- 1、Please do not hesitate to contact us if you have any questions.
- 2、Please fill in the Technical Support Request Form and email the form to support@yfxbio.com within 3months of receipt of the goods, and we will response to your inquiry within 24 hours.
- 3、The after-sale service channel will be shut down after 3 months.
- 4、Product has not been fully validated for medical applications. For research use only.