

---

## Goat Anti-Rabbit IgG (H+L), HRP

Catalog Number: YFSA02; Size: 100 $\mu$ L

---

### Background:

Yfxbio® secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA.

We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application.

---

### Product Information:

**Product name:** Goat Anti-Rabbit IgG (H+L), HRP

**Product Categories:** Secondary Antibody

**Form:** Liquid

**Tested applications:** WB, IHC-P, ELISA

**Dilution Ratio:** WB: 1:5000-1:20000; IHC: 1:500-1:1000; ELISA: 1:5000-1:20000

**Species reactivity:** Rabbit

**Host Species:** Goat

**Conjugation:** HRP

**Purity:** Affinity purification

**Immunogen:** Mouse IgG

---

### Storage:

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

---

### Services & Support:

- 1、Please do not hesitate to contact us if you have any questions.
- 2、Please fill in the Technical Support Request Form and email the form to [support@yfxbio.com](mailto:support@yfxbio.com) within 3months of receipt of the goods, and we will response to your inquiry within 24 hours.
- 3、The after-sale service channel will be shut down after 3 months.
- 4、Product has not been fully validated for medical applications. For research use only.