

Anti-beta Tubulin Rabbit Monoclonal Antibody

Catalog Number: YFRMA0053; Size: 100μL

Background:

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Product Information:

Product name: Anti-beta Tubulin Rabbit Monoclonal Antibody

Isotype: IgG

Form: Liquid

Tested applications: WB, IHC-P, ICC/IF, FC

Dilution Ratio: WB: 1:5000-1:50000; IHC: 1:500-1:2000; IF: 1:50; FC: 1:100

Species reactivity: Human, Mouse, Rat

Host Species: Rabbit

Purity: Affinity purification

Immunogen Information:

Immunogen: A synthesized peptide derived from human beta Tubulin.

Calculated molecular weight: 50kDa

Observed molecular weight: 55kDa

Swiss Prot: P07437

Gene ID : 203068

Alternative Names: Beta 4 tubulin; TBB5; TUBB2; TUBB2A; tubulin beta 2A; beta tubulin; tubulin beta chain; Tubulin beta 2b; Alpha tubulin; Alpha-tubulin ubiquitous.

Storage:

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Services & Support:

- 1、Please do not hesitate to contact us if you have any questions.
- 2、Please fill in the Technical Support Request Form and email the form to support@yfxbio.com within 3months of receipt of the goods, and we will response to your inquiry within 24 hours.
- 3、The after-sale service channel will be shut down after 3 months.
- 4、Product has not been fully validated for medical applications. For research use only.